


Initial: 2/13/08
 Reviewed/revised:
 Revision:

**MILWAUKEE COUNTY EMS
 PRACTICAL SKILL
 ECG TRANSMISSION**

Approved by: Ronald Pirrallo, MD, MHSA
 Signature:
 Page 1 of 1

Purpose:		Indications:	
To transmit 4- lead and 12-lead electrocardiograms		Any patient who requires cardiac monitoring	
Advantages:		Disadvantages:	Complications:
Transmits ECG to medical control and enables faxing 12-lead ECG to receiving hospital		None	None
		Contraindications:	None



NOTES:

- The 12-lead will remain in the Rosetta storage as long as the Rosetta device is powered on or for 20 minutes after the device turns itself off
- Paramedics can move the Rosetta and UHF radio to transmit from an area with better reception after disconnecting the Rosetta from the Zoll